#### Demonstration of hkl

Demonstrate the <a href="hkl">hkl</a> (https://people.debian.org/~picca/hkl/hkl.html) package for diffractometer computations.

First, consider you have a sample of "demo" with known unit cell parameters, mounted on a 6-circle (non kappa) diffractometer:

This example was derived from the unit test code: <a href="https://repo.or.cz/hkl.git/blob/HEAD:/tests/bindings/polarisation.py">https://repo.or.cz/hkl.git/blob/HEAD:/tests/bindings/polarisation.py</a>) with settings from <a href="https://repo.or.cz/hkl.git/blob/refs/heads/next:/tests/bindings/crystal.ini">https://repo.or.cz/hkl.git/blob/refs/heads/next:/tests/bindings/crystal.ini</a> (https://repo.or.cz/hkl.git/blob/refs/heads/next:/tests/bindings/crystal.ini)

# these other parameters from crystal.ini

Wavelength 1.62751693358

R0 0 0.0 8.0 0.0 0 1 0.0 22.31594 89.1377 0.0 0.0 45.15857 R1 1 0.0 12.0 1.0 0 1 0.0 34.96232 78.3139 0.0 0.0 71.8007

Engine hkl

Mode constant\_phi\_vertical

PsiRef not available in current engine mode

#### Before you start ...

```
In the console window:
```

```
# make sure to use the bash shell
bash

# use python3
source /APSshare/anaconda3/x86_64/bin/activate

# setup the custom hkl environment
(base) jemian@wow ~ $ . /APSshare/linux/64/hkl-5/hkl_environment.sh

# start jupyter
jupyter-notebook
```

You'll see jupyter start reporting messages in the console:

```
[I 11:21:25.136 NotebookApp] Serving notebooks from local directory: /home
/JEMIAN
[I 11:21:25.136 NotebookApp] The Jupyter Notebook is running at:
[I 11:21:25.137 NotebookApp] http://localhost:8888/?token=_____very_long_
unique_hexadecimal_code_here____
[I 11:21:25.137 NotebookApp] Use Control-C to stop this server and shut do
wn all kernels (twice to skip confirmation).
[C 11:21:25.143 NotebookApp]

To access the notebook, open this file in a browser:
    file:///run/user/970/jupyter/nbserver-24247-open.html
Or copy and paste one of these URLs:
    http://localhost:8888/?token=_____very_long_unique_hexadecimal_co
de_here_____
... more to come, ignore until you're done
```

# imports

```
In [1]: from collections import namedtuple
    from gi.repository import GLib
    import gi
        gi.require_version("Hkl", "5.0")
        from gi.repository import Hkl
        import math
        from numpy import (linspace, empty)
        from numpy.linalg import inv, norm
        import os
        import pyRestTable
        import unittest
H NUL = 12 3984244 # voltage*wavelength product angstrom * keV
```

## declare the sample lattice parameters

```
In [2]: sample = Hkl.Sample.new("demo")
        lattice = Hkl.Lattice.new(
            4.542,
                                  # a, angstrom
            16.955,
                                  # b, angstrom
            7.389,
                                  # c, angstrom
            math.radians(90.0), # alpha, radians
            math.radians(90.0), # beta, radians
            math.radians(90.0) # gamma, radians
        sample lattice set(lattice)
In [3]: wavelength = 1.62751693358
        energy = H NII / wavelength
In [4]: detector = Hkl Detector factory new(A) # TODO: what other values?
In [5]: # factories will be used to make the *hkl* objects we'll need
        factories = Hkl.factories()
        # list the diffractometers known to this version of *hkl*
        print(f"diffractometers: {sorted(factories.keys())}")
        # our diffractometer is E6C: Eulerian 6-circle (not kappa or Petra P23)
        # define a "factory" since we'll use it later
        diffractometer type = "E6C"
        e6c_factory = factories[diffractometer_type]
        # axes: MU, OMEGA, CHI, PHI, GAMMA, DELTA
        nrint(f"{diffractometer type} axes. {e6c factory create new geometry() axis names
        diffractometers: ['E4CH', 'E4CV', 'E6C', 'K4CV', 'K6C', 'PETRA3 P09 EH2', 'PETRA
        3 P23 4C', 'PETRA3 P23 6C', 'SOLEIL MARS', 'SOLEIL SIRIUS KAPPA', 'SOLEIL SIRIUS
        TURRET', 'SOLEIL SIXS MED1+2', 'SOLEIL SIXS MED2+2', 'SOLEIL SIXS MED2+3', 'SOLE
        IL SIXS MED2+3 v2', 'TwoC', 'ZAXIS']
        E6C axes: ['mu', 'omega', 'chi', 'phi', 'gamma', 'delta']
In [6]: # sample orientation and reflection 0 (0 8 0)
        angles = [0.0, 22.31594, 89.1377, 0.0, 0.0, 45.15857]
        geometry = e6c factory.create new geometry()
        geometry.axis values set(angles, Hkl.UnitEnum.USER)
        geometry.wavelength set(wavelength, Hkl.UnitEnum.USER)
        or\theta = sample add reflection(geometry detector <math>\theta + \theta = \theta
In [7]: # add reflection or1 (0 12 1)
        angles = [0.0, 34.96232, 78.3139, 0.0, 0.0, 71.8007]
        geometry.axis_values_set(angles, Hkl.UnitEnum.USER)
        orl = sample add reflection(geometry_detector_0_12_1)
```

```
In [8]: # Helper methods
                               def hkl_matrix_to_numpy(m):
                                             M = empty((3, 3))
                                             for i in range(3):
                                                          for j in range(3):
                                                                        M[i, j] = m.get(i, j)
                               def from_numpy_to_hkl_vector(v):
                                             V = Hkl.Vector()
                                             V.init(v[0], v[1], v[2])
                                             return V
  In [9]: # compute UB with or0 and or1
                                sample.compute UB busing levy(or0, or1)
                               IIB = hkl matrix to numnv(sample IIB get())
In [10]: | # compute angles for reciprocal lattice vector h, k, l
                               engine_name = "hkl"
                                engines = e6c_factory.create_new_engine_list()
                               engines.init(geometry, detector, sample)
                                engine = engines.engine_get_by_name(engine_name)
                                print(engine.modes_names_get())
                               # pick our mode
                               engine current mode set("constant nhi vertical")
                               ['bissector_vertical', 'constant_omega_vertical', 'constant_chi_vertical', 'cons
                              tant_phi_vertical', 'lifting_detector_phi', 'lifting_detector_omega', 'lifting_detector_mu', 'double_diffraction_vertical', 'bissector_horizontal', 'double_diffraction_horizontal', 'psi_constant_vertical', 'psi_constant_horizontal', 'constant_vertical', 'psi_constant_horizontal', 'constant_vertical', 'psi_constant_horizontal', 'constant_vertical', 'constant_vertic
                               nt mu horizontal']
Out[10]: 1
```

That defines our diffractometer setup.

```
In [11]: def calc(hkl):
           # assumes these are known: sample, geometry, detector, engine
           solutions = engine.pseudo_axis_values_set(hkl, Hkl.UnitEnum.USER)
           for i, s in enumerate(solutions.items()):
              values = s.geometry_get().axis_values_get(Hkl.UnitEnum.USER)
              # print(f"solution \{i+1\}: {values\}")
           first_solution = solutions.items()[0]
           values = first_solution.geometry_get().axis_values_get(Hkl.UnitEnum.USER)
           # print("picking first one")
           return values
       table = pyRestTable.Table()
       table.labels = "h k l".split()
       table.labels += geometry.axis names get()
       hkl = (0, 1, 1)
       angles = calc(hkl)
       # print(f"{hkl}: {angles}")
       table.addRow(list(hkl) + angles)
       for reflection in (or0, or1):
          hkl = reflection.hkl_get()
           angles = calc(hkl)
           # print(f"{hkl}: {angles}")
           table.addRow(list(hkl) + angles)
       nrint(tahle)
       ===
             l mu omega
       h
                                    chi
          k
                                                    phi gamma delta
       ===
              1 0.0 3.4824458166048444 22.712897698011936 0.0 0.0 13.799774663132
       0
       288
       0.0 8.0 0.0 0.0 22.31594087562736 89.13769999977886 0.0 0.0
                                                            45.158571742842
       376
       0.0 12.0 1.0 0.0 34.963469180020944 78.33265876350477 0.0 0.0
                                                           71.800704217914
       22
```

### step scan from or0 to or1

```
In [12]: num points = 5
         # get the sequence of hkl values for each Miller index
         # reference by name, not position
         # these are named tuples
         _s = or0.hkl_get()
         _f = or1.hkl_get()
         _h = linspace(_s.h, _f.h, num_points)
         _k = linspace(_s.k, _f.k, num_points)
_l = linspace(_s.l, _f.l, num_points)
         table = pyRestTable.Table()
         table.labels = "h k l".split()
         table.labels += geometry.axis names get()
         for hkl in zip(_h, _k, _l):
            angles = calc(hkl)
             #print(f"({hkl[0]:6g} {hkl[1]:6g} {hkl[2]:6g}): {angles}")
            table.addRow(list(hkl) + angles)
        nrint(tahle)
        ___ ___ ___ ___
        h k l mu omega
                                          chi
                                                            phi gamma delta
        0.0\ 8.0\ 0.0\ 0.0\ 22.31594087562736\ 89.13769999977886\ 0.0\ 0.0\ 45.158571742842
        0.0 9.0 0.25 0.0 25.15499515419808 85.49774632991891 0.0 0.0
                                                                       51.295005992443
        0.0 10.0 0.5 0.0 28.214868205373786 82.60530989054361 0.0 0.0
                                                                       57.776188360860
                                                                       64.602307648447
        0.0 11.0 0.75 0.0 31.482805809422153 80.26265823893127 0.0 0.0
        0.0 12.0 1.0 0.0 34.963469180020944 78.33265876350477 0.0 0.0
                                                                      71.800704217914
        22
```